

Sr-04-RC

**STRONTIUM-90 IN WATER CONTAINING OTHER
RADIOISOTOPES BY Cerenkov COUNTING**

Contact Person(s) : Salvatore C. Scarpitta

APPLICATION

The following procedure is used in the EML Quality Assessment Program (QAP; Sanderson and Greenlaw, 1996) for water or gross alpha/beta samples containing ^{90}Sr . This procedure was developed by Jimmy Chang, Institute of Nuclear Energy Research, Taiwan, and was tested at EML using NIST $^{90}\text{Sr}/\text{Y}$ reference materials and ^{90}Sr contained in QAP water samples. It allows for the rapid determination of ^{90}Sr ($\beta_{\text{max}} = 0.546 \text{ MeV}$) via its progeny, ^{90}Y ($\beta_{\text{max}} = 2.27 \text{ MeV}$) in aqueous solutions by Cerenkov counting (see Procedure Ba-01-R; Scarpitta and Fisenne, 1996).

Cerenkov counting is applicable for β particles with maximum energies $> 0.263 \text{ MeV}$. Alpha and gamma emitting nuclides are not detected. Immediately following separation of ^{90}Y from ^{90}Sr , a baseline count rate, C_b , is obtained to quantify any Cerenkov contribution by other nuclides that may be present in the sample.

Yield recovery is determined by adding 30 mg of Sr^{+2} carrier to the sample or by gamma counting the added ^{85}Sr tracer. The sample is Cerenkov counted at any two time intervals, t_1 and t_2 , to measure the ingrowth of ^{90}Y above C_b . Strontium-90 is calculated from the ingrowth of ^{90}Y . The Cerenkov counting efficiency for ^{90}Y in water is about 65-70%, whereas that of ^{90}Sr is $< 0.3\%$ for a 0-15 keV Cerenkov counting window. For a 20-min count time, the detection limit is about 6 mBq mL^{-1} (0.16 pCi mL^{-1}) or 0.35 dpm mL^{-1} with a relative standard deviation $< 3\%$.

SPECIAL APPARATUS

Packard Tri-Carb 2250-CA liquid scintillation spectrometer

SPECIAL REAGENTS

1. NIST traceable $^{90}\text{Sr}/\text{Y}$ reference standard of known activity (about 1000 dpm g^{-1}) to determine the ^{90}Y Cerenkov counting efficiency.
2. TRU Resin - Eichrom Industries, 8205 Cass Ave., Suite 107, Darien, IL 606651.

DETERMINATION

1. Dispense about 8-16 Bq of ^{90}Sr gravimetrically into either a 20 mL low ^{40}K borosilicate glass or plastic scintillation vial (see **Note 1**).
2. Add 10 mL of deionized water.
3. Prepare a blank using 10 mL of water.
4. Count both samples three times for 10-min each using the net average count rate to determine the ^{90}Y counting efficiency, $E^{90\text{Y}}$ (counts min^{-1} dpm $^{-1}$).
5. Preset the Packard Tri-Carb to Protocol 4 for Cerenkov determinations (Scarpitta and Fissenne, 1996). (**Note:** The Cerenkov counting window is typically 0-15 keV, although the full window, 0-2000 keV, may be used with a 50% increase in background.)
6. Obtain a sample containing an unknown amount of ^{90}Sr .
7. Add 30 mg of Sr^{++} carrier (as nitrate) to the sample for yield recovery.
8. Prepare an identical vial containing water as a Sr^{++} reference standard.
9. Reduce a premeasured amount, M_g (g), of sample to be tested to 10 mL to improve counting statistics.

Note :

1. A wavelength shifter, ANSA (7-Amino 1,3 Naphthalene di-Sulfonic Acid) can be used to enhance the Cerenkov counting efficiency but is not recommended if strontium yield recovery is to be determined gravimetrically. Nuclides that produce a Cerenkov signal in 25 mM ANSA are shown in Figure 1 (see Scarpitta and Fisenne, 1996).

SEPARATION

1. Separate the ^{90}Y from the ^{90}Sr by either oxalate precipitation (see Procedure Sr-03-RC) or EiChrom's TRU Spec extraction chromatographic resin. Record the separation date and time, t_0 .
2. Obtain a net baseline count rate for C_b immediately following ^{90}Y separation, using Protocol No. 4 on the Packard Tri-Carb 2250 CA counter and the Cerenkov counting window (0-15 keV).

CALCULATIONS

1. Recount the ^{90}Sr fraction three times a day over a 2-day period using the count rates (counts min^{-1}) CT_1 , CT_2 and CT_3 to calculate the ^{90}Sr activity in Step 2. The times t_1 , t_2 and t_3 are the number of hours after ^{90}Y separation in Step 1 of **Separation**.
2. Use the ^{90}Sr calculation as follows when ^{89}Sr is not present in the sample. (**Note:** A Basic computer program is provided in the Appendix to perform the ^{90}Sr calculations.)

$$A_1 \text{ } ^{90}\text{Sr} \text{ (Bq kg}^{-1}\text{)} = \frac{(CT_2 - CT_1) - C_b}{60 \times M_s \times E \text{ } ^{90}\text{Y} \times [\exp\{-\lambda(t_1 - t_0)\} - \exp\{-\lambda(t_2 - t_0)\}]} \quad (1)$$

$$A_2 \text{ } ^{90}\text{Sr} \text{ (Bq kg}^{-1}\text{)} = \frac{(CT_3 - CT_1) - C_b}{60 \times M_s \times E \text{ } ^{90}\text{Y} \times [\exp\{-\lambda(t_1 - t_0)\} - \exp\{-\lambda(t_3 - t_0)\}]} \quad (2)$$

where:

λ = ^{90}Y decay constant - 0.01083 h^{-1}

$E^{90}\text{Y}$ = ^{90}Y Cerenkov counting efficiency ($\text{counts min}^{-1} \text{ dpm}^{-1}$)

M_s = mass of sample (kg)

3. Obtain the average of the two ^{90}Sr activity concentrations, A_1 and A_2 from Step 2.
(**Note:** A third count could be obtained with Equation 2 if modified accordingly.)
4. Using the sample vial and the Sr^{+2} reference standard, precipitate the strontium as the carbonate, filter, dry and weigh each to obtain the yield recovery. Correct the value obtained in Step 3, dividing by the yield recovery factor. (**Note:** Gamma emitting ^{85}Sr can be added to the sample in Step 7 of **Determination** instead of Sr^{+2} .)

REFERENCES

Sanderson, C. G. and P. Greenlaw

“Semi-Annual Report of the Department of Energy, Office of Environmental Management, Quality Assessment Program”

USDOE Report EML-581, July (1996)

Scarpitta, S. C. and I. M. Fisenne

“Cerenkov Counting as a Complement to Liquid Scintillation Counting”

Appl. Radiat. Isot., 47, 795-800 (1996)

Scarpitta, S. C. and I. M. Fisenne

“Calibration of a Liquid Scintillation Counter for Alpha, Beta and Cerenkov Counting”

USDOE Report EML-583, July (1996)

Cerenkov Efficiency in 10 ml of 25 mM ANSA (Plastic Vials)




Figure 1. Cerenkov efficiency in 10 mL of 25 mM ANSA (plastic vials).

APPENDIX

Basic Computer Program

```

10  REM "SR90CER.BAS" BY J.CHANG METHOD
20  REM CALCULaTE 90SR ACTIVITY(DPM) FOR A SAMPLE USING CERENKOV COUNTING
50  TITLE$="--Test of Cerenkov Method Using 895 dpm of NIST 90-Sr Standard"
60  COUNTIME$="10 min"
100 REM ENTER TIME DATA HERE
120 DELT1=.75 ' HRS
122 DELT2=19.22
124 DELT3=111.6
130 REM ENTER COUNT RATES HERE
132 CT1=36.5
134 CT2 = 159.1
136 CT3=514.4001
150 LAMBD2 = .01083 ' HR-1
160 EFFY90=.82 ' CPM/DPM
200 REM CALCULATE SR90 ACTIVITY(DPM) HERE
205 DF1 = EXP ((-LAMBD2*DELT1)) - EXP ((-LAMBD2*DELT2))
210 A1SR90 = (CT2-CT1) / (DF1*EFFY90)
255 DF2 = EXP ((-LAMBD2*DELT1)) - EXP ((-LAMBD2*DELT3))
260 A2SR90 = (CT3-CT1)/(DF2*EFFY90)
500 REM PRINT RESULTS
510 CLS : SCREEN 2 : KEY OFF
520 PRINT DATES;" ";
525 PRINT "STRONTIUM-90 RESULTS BY Sr-SPEC and CERENKOV COUNTING" : PRINT
530 PRINT TAB(10);TITLE$ : PRINT
535 PRINT TAB (15); "Count Time = ";COUNTIME$; TAB (40) ;"Y-90 Efficiency = ";EFFY90:
    PRINT : PRINT
540 PRINT "Data" : PRINT
550 PRINT TAB (1); "Del Ti" ;TAB (10); "Cnt Ti" ;TAB (30); "Del T2" ;TAB (40); "Cnt T2" ;TA (60) ;"Del T3" ;TAB (70);
    "Cnt T3"
555 PRINT TAB (1);" (hr) ";TAB (10);" (cpm) ";TAB (30);" (hr) ";TAB (40);" (cpm) ";TAB (60); (hr) ";TAB(70);" (cpm)":
    PRINT
557 PRINT TAB(1) ;DELT1;TAB(10) ;CT1;TAB(30) ;DELT2;TAB(40) ;CT2;TAB(60) ;DELT3;TAB( 0) ;CT3
560 PRINT : PRINT : PRINT "Sr-90 Results";
570 PRINT TAB (30) ;A1SR90;" dpm" ;TAB(60) ;A2SR90;" dpm"
575 PRINT TAB(1) ;"Obs/Exp";TAB(30) ;A1SR90/895;" ";TAB(60) ;A2SR90/895;
580 PRINT : PRINT : PRINT

```

02-21-1995 STRONTIUM-90 RESULTS BY Sr-SPEC and CERENKOV COUNTING

Test of. Cerenkov Method Using 895 dpm of NIST 90-Sr Standard

Count Time -- 10 min Y-90 Efficiency -- .82

Data

Del T1 (hr)	Cnt T1 (cpm)	Del T2 (hr)	Cnt T2 (cpm)	Del T3 (hr)	Cnt T3 (cpm)
.75	36.5	19.22	159.1	111.6	514.4001
Sr-90 Results		831.421 dpm		840.6198 dpm	
Obs/Exp		.928962		.9392401	
Ok					